



The Evolution of Blood Safety Mitigation Measures

IPFA 2016, Lisbon

Adonis Stassinopoulos, Ph.D
VP Global Scientific Affairs and Research

CERUS

Outline

- Unsustainability of Testing as the main method of protection against new pathogens
- PI vs. NAT
- Introduction to INTERCEPT 10+ year experience
- Applications of PI in residual safety and new safety issues
 - Bacteria
 - Emerging Pathogens


Conclusions

Emerging Pathogens: Is testing a sustainable solution?

Development time, cost, continual addition of new pathogens

25 years of testing = protection against 9 agents

HIV, Hepatitis B, Hepatitis C, HTLV, bacteria, West Nile virus, T. Cruzi, Syphilis, Leukocytes



2003

Pathogen Inactivation

Pathogen Inactivation vs. Testing

Pathogen Inactivation



- One intervention per blood product
- Proactive
- Does not require prior knowledge of the pathogen, or additional development time
- If properly designed, ensures robust, broad inactivation
- Simpler business model

Nucleic Acid, or Ab Testing

- Developed per pathogen, and can be bundled, i.e. triplex vs. singlex
- Reactive
- Requires development and validation for sensitivity and specificity with characterized samples
- Sometimes needs to be used in addition to orthogonal assays, ie Ab and NAT
- Business model required for each pathogen

INTERCEPT Blood System

Pathogen Reduction System



Donated plasma or platelet component(s) may contain harmful agents such as bacteria, viruses, protozoans, and/or white blood cells.



When pathogens are unable to replicate, they are considered "inactivated" and cannot infect patients.




Pathogen-reduced component can then be transfused into the patient

¹INTERCEPT Blood System for Platelets [Package Insert]. Concord, CA: Cerus Corporation; 2016.

²INTERCEPT Blood System for Plasma [Package Insert]. Concord, CA: Cerus Corporation; 2015.

10+ Years of Routine Global Use

>100 Blood Centers in 25 Countries with kits sold to produce over 3,500,000 INTERCEPT platelet and plasma units




Last update: May 2016

MKT-EN 00156 v4

Select Regulatory Approvals

Status as of May 2016



Approvals shown do not represent a comprehensive list

* For first center approval

Approved Claims by Reg. Authorities

- Use of INTERCEPT platelets and plasma without patient population exclusions
- Replacement of gamma irradiation for the prevention of TA-GvHD
 - Exemptions By AABB for Accredited Labs Internationally
 - Accepted by AABB as an alternative in the US
- Replacement of CMV testing
- Extension of platelet shelf life to 7 days, when that is allowed
- Satisfies AABB requirement 5.1.5.1, as an alternative to BD
- Approved PI claims

INTERCEPT Blood Systems for All components



for Platelets



for Plasma




for RBC¹



¹The INTERCEPT System for RBC is not approved for use

Common Mechanism of Action and Target

Different compounds for PC/PL and RBC




Agents Inactivated by INTERCEPT (Amotosalen)

ROUTINELY TESTED AGENTS




ENVELOPED VIRUSES

HIV-1
HIV-2
HBV
HCV
HTVL-I
HTLV-II



SPIROCHETES

Treponema pallidum



BACTERIA

CERUS'S COMPREHENSIVE SPECTRUM PROTECTION



ENVELOPED VIRUSES

HIV-1
HIV-2
HBV
HCV
HTVL-I
HTLV-II

BVDV
CMV
WNV
SARS
Vaccinia¹
Chikungunya
Dengue²
Influenza A

DH BV



NON-ENVELOPED VIRUSES

Bluetongue virus, type 11
Simian Adenovirus-15
Feline calicivirus
Parvovirus B19
Human adenovirus 5




GRAM-NEGATIVE BACTERIA

Klebsiella pneumoniae
Yersinia enterocolitica
Escherichia coli
Pseudomonas aeruginosa
Salmonella choleraesuis
Enterobacter cloacae
Serratia marcescens
Anaplasma phagocytophilum
*Orientia tsutsugamushi*³



GRAM-POSITIVE BACTERIA

Staphylococcus epidermidis
Staphylococcus aureus
Streptococcus pyogenes
Listeria monocytogenes
Corynebacterium minutissimum
Bacillus cereus (vegetative)
Lactobacillus sp.
Bifidobacterium adolescentis
Propionibacterium acnes
Clostridium perfringens



SPIROCHETES

Treponema pallidum
Borrelia burgdorferi



PROTOZOA

Trypanosoma cruzi
Plasmodium falciparum
Leishmania sp.
Babesia microti



LEUKOCYTES

T-cells

(1) Sampson-Johannes A, et al. 2003. *Transfusion*. 43:83A; (2) Lam S, et al. *Transfusion* 2007;47:134A; (3) Rentas F. *Transfusion* 2004;44:104A.


Agents Inactivated by INTERCEPT (Amustaline)

ROUTINELY TESTED AGENTS




ENVELOPED VIRUSES

HIV-1
HIV-2
HBV
HCV
HTVL-I
HTLV-II



SPIROCHETES

Treponema pallidum



BACTERIA

CERUS'S COMPREHENSIVE SPECTRUM PROTECTION



ENVELOPED VIRUSES

HIV-1
HIV-2
HBV
HCV
HTVL-I
HTLV-II

BVDV
CMV
WNV
SARS
Vaccinia¹
Chikungunya
Dengue²
Influenza A

DH BV



NON-ENVELOPED VIRUSES

Bluetongue virus, type 11
Simian adenovirus 15
Feline calicivirus
Parvovirus B19
Human Adenovirus 5




GRAM-NEGATIVE BACTERIA

Klebsiella pneumoniae
Yersinia enterocolitica
Escherichia coli
Pseudomonas aeruginosa
Salmonella choleraesuis
Enterobacter cloacae
Serratia marcescens
Anaplasma phagocytophilum
*Orientia tsutsugamushi*³



GRAM-POSITIVE BACTERIA

Staphylococcus epidermidis
Staphylococcus aureus
Streptococcus pyogenes
Listeria monocytogenes
Corynebacterium minutissimum
Bacillus cereus (vegetative)
Lactobacillus sp.
Bifidobacterium adolescentis
Propionibacterium acnes
Clostridium perfringens



SPIROCHETES

Treponema pallidum
Borrelia burgdorferi



PROTOZOA

Trypanosoma cruzi
Plasmodium falciparum
Leishmania sp.
Babesia microti



LEUKOCYTES

T-cells

(1) Sampson-Johannes A, et al. 2003. *Transfusion*. 43:83A; (2) Lam S, et al. *Transfusion* 2007;47:134A; (3) Rentas F. *Transfusion* 2004;44:104A.



Hemovigilance Programs

Demonstrated safety in routine use

**>300,000 INTERCEPT Platelets
Evaluated in Routine Use**

**French National
Hemovigilance¹**

**Swiss National
Hemovigilance²**

**Multicenter
Cerus HV³**

**# INTERCEPT Platelet
Transfusions**

180,762

130,843

19,175

**# Patients Receiving
INTERCEPT Platelets**

~30,000

~20,000

4,067

INTERCEPT ATR Rate

~0.3%⁴

~0.3%⁵

~0.6%

Conventional ATR Rate

~0.3%⁴

~0.4%⁵

NA


(1) French National Agency for Medicine and Health Product Safety/ANSM, Hemovigilance Activity Reports, 2006 - 2014. (2) SwissMedic Haemovigilance Annual Reports, 2010 - 2014. (3) Knutson F et al. Vox Sanguinis 2015. (4) French National Agency for Medicine and Health Product Safety/ANSM, Hemovigilance Activity Reports, 2009-2011. (5) SwissMedic Haemovigilance Annual Report 2014; conventional : 2008-2010, INTERCEPT: 2011-2014.

BACTERIA

FDA Revised Draft Guidance for Bacterial Safety¹

Recognizes bacterial risk in platelets, includes pathogen reduction as option to help mitigate risk

- Addresses industry concerns on the high risk of bacterial contamination in platelets.
- Includes pathogen reduction as an alternative to testing, strengthens language for day 4 and 5 secondary testing requirements.
- Draft Guidance includes:
 - Blood centers - must either:
 - Perform **pathogen reduction**, or
 - Perform primary bacterial testing (use of culture detection method on a sample taken ≥24 hours post-collection)
 - Hospital transfusion service – can use either:
 - **Pathogen reduced** platelets – no further measures necessary
 - Conventional platelets, already negative in primary bacterial testing -- must perform secondary testing at day 4 and 5 per methods in guidance document



1. "Bacterial Risk Control Strategies for Blood Collection Establishments and Transfusion Services to Enhance the Safety and Availability of Platelets for Transfusion," FDA Draft Guidance for Industry, March 2016.

High Log Inactivation for Broad Spectrum of Bacteria in Platelet Concentrates

Gram-positive (Aerobes and Anaerobes)	PC/PAS	PC/PL	Gram-negative	PC/PAS	PC/PL
Staphylococcus epidermidis	> 6.6	>7.4 ¹	Escherichia coli	> 6.4	≥7.3
Staphylococcus aureus	6.6	>7.6	Serratia marcescens	> 6.7	>7.1
Streptococcus pyogenes	> 6.8	>6.1	Klebsiella pneumoniae	> 5.6	≥6.7
Listeria monocytogenes	> 6.3	>6.6	Pseudomonas aeruginosa	4.5	≥6.8
Corynebacterium minutissimum	> 6.3	>6.4	Salmonella choleraesuis	> 6.2	>5.9
Bacillus cereus (vegetative)	> 6.0	≥5.6	Yersinia enterocolitica	> 5.9	>7.3
Lactobacillus sp	> 6.9	>6.1	Enterobacter cloacae	5.9	≥6.0
Bifidobacterium adolescentis	> 6.5	-	Spirochetes		
Propionibacterium acnes	> 6.7	>6.7	Treponema pallidum	6.8 - 7.0	-
Clostridium perfringens	> 7.0	>6.0	Borrelia burgdorferi	> 6.8	-

- **Robust Inactivation of Bacteria is critical to protect against them, as bacteria can grow in blood products**

- ¹ In blue latest claims that have not been reviewed by TUV yet

No Sepsis Fatalities with INTERCEPT Platelets when Used in Routine

French and Swiss National hemovigilance data

	Year	Conventional Platelets		INTERCEPT Platelets	
		Units Transfused (n)	Transfusion Transmitted Infections (Fatalities)	Units Transfused (n)	Transfusion Transmitted Infections (Fatalities)
 FRANCE	2006	231,853	4 (0)	6,420	0 (0)
	2007	232,708	9 (2)	15,393	0 (0)
	2008	239,349	6 (1)	15,544	0 (0)
	2009	241,634	9 (0)	21,767	0 (0)
	2010	253,149	2 (1)	22,632	0 (0)
	2011	267,785	3 (1)	22,392	0 (0)
	2012	275,834	7 (2)	24,849	0 (0)
	2013	278,234	4 (1)	25,089	0 (0)
	2014	278,788	5 (0)	26,676	0 (0)
 SWITZERLAND	2010	29,900	1 (0)	0	0 (0)
	2011	6,600	0	26,500	0 (0)
	2012	0	0	34,265	0 (0)
	2013	0	0	34,750	0 (0)
	2014	0	0	35,328	0 (0)
Total		2,335,834	50 (8)	311,605	0 (0)

French HV data: 2006-2011: Sweeney J, Lozano M. Platelet Transfusion Therapy. Bethesda: AABB Press, 2013.

2012-2014: French National Agency for Medicine and Health Product Safety/ANSM, Hemovigilance Activity Reports.

Swiss HV data: SwissMedic Haemovigilance Annual Reports, 2010-2014.

EMERGING PATHOGENS

Emerging Pathogens

The risk of spreading pathogens – Chikungunya, Dengue, Zika

Chikungunya Cases Spike In New York And New Jersey

Reuters
Posted: 07/30/2014 1:13 pm EDT | Updated: 09/29/2014 6:59 am EDT



Hawaii mayor declares state of emergency to address dengue fever outbreak



Zika virus is an illness transmitted to people through bites from mosquitoes that spread dengue and chikungunya viruses.

Chagas and chikungunya in Texas

Chagas and Chagas disease in those areas may be a risk in spreading fatal

2014

Hawaii Department of Health Reports First Case Of Microcephaly Caused By Zika Virus In The U.S. (Updated)

Andrew Lipnik
1/10/16 2:22pm - Filed for ZIKA VIRUS



Just days after the Center for...
warnings to...

Published On: Thu, Sep 12th, 2013

Hometown News / Outbreak News | By Robert Herriman

Florida Blood Bank, OneBlood, Temporarily Suspends Blood Donations From Martin And St. Lucie Counties Due To Dengue Fever

Emerging Pathogens

Current mitigations: Pathogen Reduction

FDA Guidance, February 2016

In areas with active Zika virus transmission

(<http://www.cdc.gov/zika/geo/index.html>), the FDA recommends that Whole Blood and blood components for transfusion be obtained from areas of the U.S. without active transmission. **Blood establishments may continue collecting and preparing platelets and plasma if an FDA-approved, pathogen-reduction device is used.** The guidance also recommends blood establishments update donor education materials with information about Zika virus signs and symptoms and ask potentially affected donors to refrain from giving blood.

Recommendations for Donor Screening, Deferral, and Product Management to Reduce the Risk of Transfusion- Transmission of Zika Virus

Guidance for Industry

This guidance is for immediate implementation.

FDA is issuing this guidance for immediate implementation in accordance with 21 CFR 10.115(g)(2) without initially seeking prior comment because the agency has determined that prior public participation is not feasible or appropriate.

Submit one set of either electronic or written comments to the Division of Dockets Management, 630 Fishers Lane, Rm. 1061, Rockville, MD 20850. Include the docket number listed in the notice of this guidance. FDA will review any comments we receive.

For more information, contact the Office of Communication, Outreach and Development, 10155 Belts Road, Bldg. 71, Rm. 3128, Silver Spring, MD 20910, or call 1-800-368-1099 or email ocod@fda.hhs.gov, or visit www.fda.gov/ocod.

www.fda.gov/ocod

For more information, contact OCOD at the phone numbers or email address listed above.

U.S. Department of Health and Human Services
Food and Drug Administration
Center for Biologics Evaluation and Research
February 2016

Current mitigations: Pathogen Reduction

WHO Guidance, February 2016

c. Pathogen reduction of blood components

Pathogen reduction technology (PRT) may be implemented. PRT is currently available for plasma and platelets, but not for whole blood or red blood cells. Different PRTs have been shown to be effective against other flaviviruses (e.g. West Nile, dengue) [10,11,20] and, in the absence of Zika-specific information, are presumed as equally effective against Zika virus.

Maintaining a safe and adequate blood supply during Zika virus outbreaks

Interim guidance
February 2016
WHO/ZIKV/HS/16.1



1. Introduction

1.1 Background

These guidelines have been developed in recognition that infection with Zika virus may present a risk to blood safety, and in consideration of the declaration on 1 February 2016 by the WHO Director-General of a Public Health Emergency of International Concern with regard to clusters of microcephaly and other neurological disorders, potentially associated with Zika virus. Currently there is

other autoimmune neurological complications was suspected during a 2013–2014 outbreak in French Polynesia and remains under investigation [6, 16].

During the Zika virus outbreak in French Polynesia between November 2013 and February 2014, a total of 1,505 healthy blood donors were tested by nucleic acid amplification technology (NAT)-based assays, with 42 (2.8%) confirmed positive for Zika virus RNA. Blood donors positive for Zika virus RNA were contacted retrospectively to investigate the occurrence of 'Zika fever'.

conjunctivitis and/or arthralgia) If the 42 donors that tested negative had a Zika fever-like fever they gave blood. No enough transfusion was [3]. However, transmission (of Zika virus) by blood (and West Nile viruses) by documented [2, 18, 22].

use by national health
services, to provide a
lines applicable to their own
text may be developed.

adequate blood
with ongoing Zika virus

through reinforcing blood
areas

virus outbreak should ideally be
ed collections in non-affected

clinical symptoms of Zika virus infection, and probable cases of sexual transmission have been described [7, 8, 14, 15].

A link between Zika virus infection during pregnancy and microcephaly in neonates is suspected and currently being investigated for causal association [12, 21]. An association of Zika virus with Guillain-Barré syndrome (GBS) and

areas. In non-affected areas, consideration may be given to defer donors who have recently visited areas with ongoing transmission of Zika virus for a period of 28 days after their departure from the area (twice the assumed maximum incubation period [9]). It is crucial that public health authorities work with the blood transfusion service (BTS) to establish mechanisms to access regular, up-to-date

The INTERCEPT Blood System

Maintaining a safe and adequate blood supply during arboviral outbreaks

INTERCEPT has been called upon to help sustain platelet availability during arboviral outbreaks:



- 2005 CHIKV epidemic on Ile de La Réunion, France¹
 - All PC were treated with INTERCEPT allowing collections of PC in the island
- CHIKV outbreaks in Guadeloupe and Martinique, French Polynesia²
 - PC were treated with INTERCEPT
 - Other components were either quarantined with enhanced PDI collection, or tested with NAT
 - CHIKV positive donations was transfused to 10 recipients without clinical symptoms³
- 2014 CHIKV and DENV outbreaks in the Caribbean region⁴
 - INTERCEPT PC components were deployed in Puerto Rico in response to the epidemic
 - Authorities excluded PI treated components from the quarantine imposed to other components
- 2016 WHO⁵ and US FDA⁶ issued guidances recommending pathogen reduction as an option to mitigate risks related to ZIKV outbreaks; AABB⁷ issued guidance as an option to mitigate CHIKV, DENV, ZIKV

1. Rasongles P, et al. Transfusion 2009;49:1083-91. 2. Musso D, et al. Transfusion. 2014 Nov;54(11):2924-30. 3. Gallian et al. Transfusion, sup 3 17A, 2015 4. Rico, S et al. Treatment Use Study of INTERCEPT Platelet Components in Response to the Chikungunya and Dengue Epidemic in Puerto Rico: TRUE Study. presented at BMT Tandem Meeting; February 2016. 5. WHO Interim guidance WHO/ZIKV/HS/16.1 February 2016. 6. FDA Zika Guidance for Industry February 2016. 7. AABB Association Bulletin #16-04. March 1, 2016.

Arbovirus Comparison

Characteristic	Dengue	Chikungunya	Zika	West Nile
Virus family	Flaviviridae	Togaviridae	Flaviviridae	Flaviviridae
Virus genus	Flavivirus	Alphavirus	Flavivirus	Flavivirus
Serotypes	4	1	1	1
RNA size (kb)	10.7	11.8	10.8	10.3
Vector	<i>Aedes aegypti</i> <i>Aedes albopictus</i>	<i>Aedes aegypti</i> <i>Aedes albopictus</i>	<i>Aedes aegypti</i> <i>Aedes albopictus</i>	<i>Culex Pipiens</i> <i>C. Tarsalis</i> <i>C. quinequefasciaatus</i>
Symptoms / illness	Acute febrile illness	Acute febrile illness	Acute febrile illness	Acute febrile illness, rash
Symptomatic: Asymptomatic	25:75	75:25	20:80 Sexually Transmitted	20:80
Illness outcome	Severe dengue – plasma leakage	Rarely severe, Arthralgias	Rarely severe, Guillian-Barré Microcephaly	Neuroinvasive Disease Meningitis Encephalitis Polyomyelitis

Pathogen Inactivation of Arboviruses in Plasma by the INTERCEPT Blood System

Arbovirus	Genus	INTERCEPT Log ₁₀ Reduction
West Nile Virus	Flavivirus	>6.8
Dengue	Flavivirus	>5.6 ¹
Chikungunya	Alphavirus	>6.4
Zika	Flavivirus	>6.6 ²

¹ Musso D. et al, Transfusion. 2014 54, 2924-30

² Aubry M. *et al.* Transfusion, 2016

Published Zika PI in INTERCEPT Plasma

TABLE 1. Detection of replicative ZIKV and ZIKV titration (log TCID₅₀/mL) in plasma samples before and after inactivation

Plasma samples	Initial viral titers	Replicative ZIKV after					Log reduction
		First passage	Second passage	Third passage	Fourth passage	Fifth passage	
A							
Preinactivation sample	6.46	+*	+	+	+	+	>6.46
Inactivated sample	-†	-	-	-	-	-	
B							
Preinactivation sample	6.63	+	+	+	+	+	>6.63
Inactivated sample	-	-	-	-	-	-	
C							
Preinactivation sample	6.61	+	+	+	+	+	>6.61
Inactivated sample	-	-	-	-	-	-	
D (control)							
Preinactivation sample	6.14	+	+	+	+	+	
Noninactivated sample	6.28	+	+	+	+	+	

* Positive immunofluorescence.

† Negative immunofluorescence.

**>6.57
Infectivity**

TABLE 2. ZIKV RNA quantitation (log copies/mL) in plasma samples before and after inactivation

Plasma sample	Initial RNA loads	RNA loads after				
		First passage	Second passage	Third passage	Fourth passage	Fifth passage
A						
Preinactivation sample	10.22	10.31	10.08	10.36	10.17	10.48
Inactivated sample	9.67	4.07	ND*	ND	ND	ND
B						
Preinactivation sample	10.41	10.30	10.17	10.26	10.22	10.42
Inactivated sample	9.44	3.63	ND	ND	ND	ND
C						
Preinactivation sample	10.11	10.43	10.19	10.30	10.20	10.65
Inactivated sample	9.41	3.87	ND	ND	ND	ND
D (control)						
Preinactivation sample	10.36	10.26	10.02	10.32	10.02	10.28
Noninactivated sample	9.90	9.89	9.99	9.95	10.18	10.01

* ZIKV RNA not detected.

**10.25
PCR**

Translating Infectivity to Detection Language¹




Zika Infectivity PI (cfu/mL) $>6.57 \log_{10}$

Zika PCR Reduction (gEq/mL) $10.25 \log_{10}$

$$\log_{10} (\text{gEq/mL}) \geq 3.68 + \log_{10} (\text{CFU/mL})$$

Inactivation in
INTERCEPT
Plasma

Donor Viremia



¹Aubry *et al.* Transfusion 2015

Zika Inactivation using the INTERCEPT System for Platelets

Conditions	N	Log Titers (pfu/mL)		Log Reduction per mL	Log Reduction
		Pre-UVA	Post-UVA		
PC in 100% PLASMA	1	4.3	< 0.0	> 4.3	> 4.3
PC in 65% PAS	2	4.0	<-1.0	>4.0	>4.9
PC in 65% PAS*	2	6.2	<-0.6	>6.2	>6.8

*Small volumes

Zika Inactivation using the INTERCEPT System for Red Blood Cells with S-303

Replicate	Log Titers (pfu/mL)			
	Control	Test	Log Reduction per mL	Log Reduction
1 LR	4.3	<-0.7	>4.3	>5.0
2 LR	4.4	<-0.7	>4.4	>5.1
3 NLR	4.3	<-0.7	>4.3	>5.0
Mean	4.3	<-0.7	>4.3	>5.0

Conducted at nominal volumes (~280 mL AS-5 RBC)

Conclusions

- Blood Safety is ensured through the use of multiple redundant measures, with testing currently being a key contributor to enhanced blood safety
- There are areas where testing appears inadequate, or cannot be developed and applied fast enough
- PI is increasingly recognized as a technology that can address such areas, as examples from the bacterial contamination and emerging pathogens have demonstrated, in routine use.
- The introduction of PI for RBC will expand the applicability of PI in enhancing blood safety

Thank you!

CERUS